

| ASTM and IP Test Methods | | ASTM D130 | ASTM D665 | ASTM D1838 | ASTM D4048 | ASTM D6138 | ASTM D7667 | ASTM D7671 | IP 112 | IP 154 | IP 227 |
|--------------------------|--|-----------|-----------|------------|------------|------------|------------|------------|--------|--------|--------|
| Accessories | | | | | | | | | | | |
| 11241-0 | Silicone carbide paper 150 grit (pack of 50) | | | | | | | | | | ● |
| 11420-0 | Support assembly (pair) | | | | | | | | | | |
| 11430-0 | Glassware silver corrosion | | | | | | | | | | ● |
| 11440-0 | Silver test strip (pack of 5) | | | | | | | | | | ● |
| 11450-0 | ASTM silver strip corrosion standard | | | | | | | | | | ● |
| 11460-0 | Silicone carbide paper 240 grit (pack of 50) | | | | | | | | | | |
| 11465-0 | Silicone carbide paper 220 grit (pack of 50) | | | | | | | | | | |
| 11470-0 | Silicone carbide paper 100 grit (pack of 50) | | | | | | | | | | |
| 11480-0 | Silicone carbide grains 150 grit 500g (pack of 5) | | | | | | | | | | ● |
| 11510-0 | Silver corrosion kit | | | | | | | | | | ● |
| 11512-0 | Copper corrosion kit for LPG | | | | | | ● | | | | |
| 11514-0 | Tongs | | | | | | | | | | ● |
| 11515-0 | Silver corrosion test kit for gasoline | | | | | | ● | | | | |
| 11516-0 | Thin silver strip corrosion kit | | | | | | | | | | |
| 11517-0 | Copper corrosion kit for gasoline and jet fuel | ● | | | | | | | | ● | |
| 11518-0 | Copper corrosion kit for diesel, fuel oil and lubricants | ● | | | | | | | | ● | |
| 11519-0 | Copper corrosion kit for grease | | | | ●* | | | | ●** | | |
| 11530-0 | Multi-strip vice | | | | | | | | | | ● |
| 11220-0 | Test specimen/holder/cover (pack of 5) | | ● | | | | | | | | |

* replace test tubes with 11000-002 cloud and pour jars, also include a vented cork for these jars

** replace test tubes with 100ml tall form beakers and remove 3 way holder and vented corks

| ASTM and IP Test Methods | | ASTM D130 | ASTM D665 | ASTM D1838 | ASTM D4048 | ASTM D6138 | ASTM D7667 | ASTM D7671 | IP 112 | IP 154 | IP 227 |
|--------------------------|---|-----------|-----------------------|------------|------------|------------|------------|------------|--------|--------|--------|
| 11240-0 | Aluminium oxide cloth 150 grit (pack of 50) | | ● | | | | | | | | |
| 11250-0 | Aluminium oxide cloth 240 grit (pack of 50) | | ● | | | | | | | | |
| 11260-0 | Baffled beaker (NACE) (pack of 4) | | ● | | | | | | | | |
| 11275-4 | Polishing and grinding rig | | ● | | | | | | | | |
| 11280-0 | Chuck body and locknut | | ● | | | | | | | | |
| 11225-0 | Synthetic seawater 30ml (pack of 10) | | ● Procedure B only | | | | | | | | |
| 11227-0 | Synthetic seawater 300ml | | ● Procedure B only | | | | | | | | |
| 11325-0 | Beaker 25ml (pack of 10) | | ● Procedure B only | | | | | | | | |
| 19301-0 | EMCOR starter kit | | | | | ● | | | | | |
| 94635-0 | Bath oil 80-120°C 20 litres | ● | ● 2 required | | ● | ● | | | | ● | |
| ASTMS12CX | Thermometer -20 to +102°C | ● | | ● | ● | ● | | ● | | ● | ● |
| ASTMST9C | Thermometer -5 to +110°C | ● | | ● | ● | ● | | ● | | ● | ● |
| Spares | | | | | | | | | | | |
| 11516-001 | Silver strip centering device (pack of 10) | | | | | | ● | | | | |
| 11516-002 | Scouring pad (pack of 10) | | | | | | ● | | | | |
| 11516-003 | Polished silver strips (pack of 20) | | | | | | ● | | | | |
| 11516-004 | Aluminium oxide paper 150 grit | | | | | | ● | | | | |
| 11240-0 | Aluminium oxide cloth 150 grit (pack of 50) | | ● | | | | | | | | |
| 11250-0 | Aluminium oxide cloth 240 grit (pack of 50) | | ● | | | | | | | | |

| ASTM and IP Test Methods | | ASTM D130 | ASTM D665 | ASTM D1838 | ASTM D4048 | ASTM D6138 | ASTM D7667 | ASTM D7671 | IP 112 | IP 154 | IP 227 |
|--------------------------|---------------------------------------|-----------|-----------------------|------------|------------|------------|------------|------------|--------|--------|--------|
| 11260-0 | Baffled beaker (NACE) (pack of 4) | | ● | | | | | | | | |
| 11325-0 | Beaker 25ml (pack of 10) | | ● Procedure B only | | | | | | | | |
| 11330-0 | Pair of test specimens | | ● | | | | | | | | |
| 11335-0 | Beaker cover and test specimen holder | | ● | | | | | | | | |
| 11200-005 | Berzelius beaker 400ml (pack of 10) | | ● | | | | | | | | |
| 11200-503 | Stirrer (pack of 2) | | ● | | | | | | | | |
| 11225-0 | Synthetic seawater 30ml (pack of 10) | | ● Procedure B only | | | | | | | | |
| 11227-0 | Synthetic seawater 300ml | | ● Procedure B only | | | ● | | | | | |
| 19300-001 | EMCOR v ring seal (pack of 16) | | | | | ● | | | | | |
| 19300-002 | EMCOR mounting nut (pack of 8) | | | | | ● | | | | | |
| 19300-003 | EMCOR sleeve (pack of 8) | | | | | ● | | | | | |
| 19300-004 | EMCOR test bearing | | | | | ● | | | | | |
| 11440-0 | Silver test strip (pack of 5) | | | | | | | ● | | | |
| 11490-0 | Silver strip glass cradle (pack of 3) | | | | | | | ● | | | |
| 11520-0 | Test tube ejector spring | | | | | | | ● | | | |
| 11550-0 | Copper test strip (pack of 30) | ● | | | | | | | ● | ● | |
| 11590-0 | Test tube (pack of 10) | ● | | | | | | | ● | ● | |
| 22190-0 | Copper test strip (pack of 30) | | | ● | ● | | | | | | |
| 11000-002 | Graduated jar and cork (pack of 10) | | | | ● | | | | | | |